
Isotype Control Mouse IgM eFluor® 450

Catalog Number: 8048-4752

Analyte Specific Reagent (ASR): Analytical and performance characteristics are not established.

Product Information

Contents: Isotype Control Mouse IgM eFluor® 450

 **Catalog Number:** 8048-4752

Clone: 11E10

Concentration: 0.2 mg/ml

Host/Isotype: Mouse IgM

Formulation: aqueous buffer, 0.09% sodium azide, may contain carrier protein/stabilizer



Temperature Limitation: Store at 2-8°C. Do not freeze. Light-sensitive material.



Batch Code: Refer to vial



Analyte Specific Reagent



Manufacturer: eBioscience, Inc., 10255 Science Center Drive, San Diego, CA 92121, USA



Use By: Refer to vial

Caution, contains Azide

Description

The immunoglobulin was generated by immunization with lipopolysaccharide (LPS) and is suitable for use as a mouse IgM isotype control.